

Cell Isolation from Nasopharyngeal Swabs for scRNA-seq - Illustrated Protocol

Introduction

This protocol can be used to isolate cells from frozen nasopharyngeal swabs. A typical good yield is 50-100,000 cells for each swab.

<https://www.protocols.io/view/human-nasopharyngeal-swab-processing-for-viable-si-5jyl8myz9g2w/v1/materials>

Samples have been collected using a nasal swab and frozen using a slow-cooling device in 90% FBS + 10% DMSO in a 1.7 mL cryovial in -80°C, stored in LN2.

Materials

Materials & Reagents for 1 swab (with overages):

- 2.5 mL of RPMI + 10 mM dithiothreitol (DTT) (made fresh)
- 3.5 mL of Accutase
- 6 mL RPMI
- 8 mL quenching buffer containing RPMI + 10% fetal bovine serum (FBS) + 4 mM EDTA
- 2 mL RPMI + 10% FBS
- 15 mL conical labeled **Tube B** containing 2 mL RPMI (previously 5 mL)
- 1.5 mL tube labeled **Tube B**, empty
- 1.5 mL tube labeled **Tube C**, with 1 mL RPMI + 10 mM DTT
- 1.5 mL tube labeled **Tube D**, with 1 mL Accutase
- 50 mL conical
- 70 µm strainer that fits 50 mL conical
- 1.5 mL tube labeled **Tube E**
- 1.5 mL tube for counting
- 1.5 mL tube for pooling cells (if doing 3 samples, need 1 tube)
- 10 µL trypan
- 1 mL PBS+ 1% BSA
- 15 mL conical for collecting viral supernatant labeled **Tube F**
- 3 cryovials for viral supernatant lysate
- 3 96-well PCR plates or 3 cryovials for bulk RNAseq lysate
- RLT + 1% 2-mercaptoethanol (BME)

Equipment


- Forceps and scissors
- Thermomixer set to 37°C, agitating at 300 rpm
- Hemocytometer slides

Set-Up

1. Record sample characteristics in Table 1

| | Study ID | Case/Control | Virus | Date of Enrollment | Sex | Age | Sample ID |
|---|----------|--------------|-------|--------------------|-----|-----|-----------|
| 1 | | | | | | | |
| 2 | | | | | | | |
| 3 | | | | | | | |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | | | | | | | |


2. Prepare 50 mL conicals with necessary reagents:
- o RPMI + 10 mM DTT
 - o Accutase
 - o RPMI + 10% FBS + 4 mM EDTA (quenching media)
 - o RPMI + 10% FBS




3. Label 2 sets of 15 mL conicals (**Tube B & Tube F**). Pre-fill **Tube B** with 2 mL RPMI




4. Label 2 sets of 96 well plates (**Tube C, Tube D, Tube E**) Pre-fill **Tube C**



5. Label 1 set of 50 mL tubes for filtering at later step (can set aside for beginning of protocol)



Complete Set-Up



Timing

Fill out timing chart at indicated steps of protocol:

Start time:

Time after preparation of single cell suspension:

Time after counting:

Time 10X Controller run started:

Time 2nd 10X controller run started:

Time Lysis Buffer added:

Procedure

Tube A (the cryovial)

2. Rapidly thaw cryovial in thermal block set to 37°C.
3. Remove swab from cryovial using forceps.
4. Place swab in **Tube B** (15 mL conical), dip briefly to rinse swab
5. Move swab from **Tube B** (15 mL conical) to **Tube C**. Trim swab handle using scissors if necessary. (continue below **Tube C**)
6. Transfer liquid in **Tube A** to **Tube B** (15 mL conical)
7. Using ~ 1 mL RPMI from **Tube B** (15 mL conical), wash **Tube A**. Collect washing in **Tube B** (15 mL conical) (continue below **Tube B**)
8. Discard **Tube A**

Tube B

9. Centrifuge 15 mL conical **Tube B** at 400g for 5 minutes at 4°C
10. Remove supernatant with serological or P1000
11. Place supernatant in **Tube F**
12. Resuspend pellet in 1 mL RPMI + 10 mM DTT

13. Transfer suspended cells from **Tube B** (15 mL) to **Tube B** (1.5 mL). Discard empty 15 mL conical.
14. Place **Tube B** 1.5 mL eppendorf on thermomixer (37°C, 300 rpm)
15. Incubate for 15 minutes
16. Centrifuge **Tube B** (1.5 mL) at 400g for 5 minutes at 4°C
17. Remove supernatant with P1000 pipette
18. Resuspend pellet in 1 mL Accutase
19. Place **Tube B** (1.5 mL) on thermomixer (37°C, 300 rpm)
20. Incubate for 30 minutes

Tube C

21. Place **Tube C** on thermomixer (37°C, 300 rpm)
22. Incubate for 15 minutes
23. Place swab in **Tube D** (continue below **Tube D**)
24. Centrifuge remaining liquid at 400g for 5 minutes at 4°C
25. Remove supernatant with P1000 pipette
26. Place supernatant in **Tube F**
27. Resuspend pellet in 1 mL Accutase
28. Place **Tube C** on thermomixer (37°C, 300 rpm)
29. Incubate for 30 minutes


Tube D

30. Place **Tube D** on thermomixer (37°C, 300 rpm)
31. Incubate for 30 minutes
32. Take 10X reagents out at this step so that they have time to equilibrate to the appropriate temperature

After **Tube B**, **Tube C**, and **Tube D** have finished their 30 minute incubation:

Note: in practice, we wait until all tubes have finished and synched up and leave tubes on incubation for longer than 30 minutes (maximum 50 minutes).

34. Place 70 μ m filter in 50 mL conical
35. Wet filter with 3 mL quenching buffer (RPMI + 10% FBS + 4 mM EDTA)
36. Pipette contents of **Tube B**, **Tube C**, and **Tube D** onto filter (Do not discard original tubes)
37. Use 1 mL of fresh quenching buffer to wash each **Tube B**, **Tube C**, and **Tube D**. Manually agitate the swab in **Tube D** in the quenching buffer with forceps to ensure full rinse
38. Add quenching buffer from washes to filter. Discard **Tubes B**, **C**, and **D**.
39. Wash filter with additional 2 mL quenching buffer
40. Discard filter, cap 50 mL conical
41. Centrifuge 50 mL conical at 400g for 10 minutes at 4°C



It will likely be very challenging to see a pellet in the 50 mL conical at this point.

42. Remove supernatant with serological. Leave ~500 μ L in the bottom of the tube

Remove supernatant carefully at this step! You don't want to pipette up your cells but you will need to add the volume you leave in this 50 mL conical and 1 mL RPMI + 10% FBS to a 1.5 mL tube in the next steps, so try not to exceed 500 μ L residual volume.



Example residual volume

43. Add 500 μ L RPMI + 10% FBS to the tube to resuspend cells.
44. Transfer resuspended cells (\sim 1 mL) to **Tube E** (1.5 mL tube)
45. Wash 50 mL conical with additional 500 μ L RPMI + 10% FBS. Transfer washing to **Tube E**
46. Centrifuge **Tube E** at 400g for 5 minutes at 4°C
47. Remove supernatant with P1000 pipette
48. Resuspend pellet in 1 mL RPMI + 10% FBS
- At this point, you should be able to see a reasonably-sized pellet!*
49. Centrifuge **Tube E** at 400g for 5 minutes at 4°C
50. Resuspend pellet in 200 μ L RPMI + 10% FBS




sized

Example pellet in Tube E after centrifuging.

Count cells from **Tube E**

51. In 1.5 mL tube, add 10 μ L trypan
52. Add 10 μ L cells from **Tube E** to 1.5 mL tube containing trypan blue.
53. Pipette to mix cells in trypan blue, transfer 10 μ L to hemocytometer port
54. Count viable cells across 4 quadrants
55. Record total cell number and cell concentration
56. **Take photo of cells at 4x**



Representative hemocytometer images at 4X

Prepare cells for 10X

57. Prepare 1.5 mL tubes for each pool
58. For each sample, add "Volume Necessary for Desired Number of Cells" from Table 3 to correct pool
59. Pipette cells to mix
60. If volume of pool is $> 100\ \mu\text{L}$, centrifuge cells at 400g for 5 minutes, then resuspend in $500\ \mu\text{L}$ PBS + 1% BSA and proceed to next step. If volume of pool is $< 100\ \mu\text{L}$, add $500\ \mu\text{L}$ PBS + 1% BSA directly to pool and proceed to next step.
61. Centrifuge cells at 400g for 5 minutes
62. Resuspend cells in $43.3\ \mu\text{L}$ PBS + 1% BSA (Total volume of cell suspension + water on page 27 of 10X Protocol)
63. Add $43.3\ \mu\text{L}$ of cell suspension directly to $31.9\ \mu\text{L}$ master mix (Step 1.2b of 10X Protocol)
64. Proceed with instructions in 10X Protocol to load the chip and run the controller (through Step 1.3)
65. **At Step 1.4f, take a picture of the GEMs in the pipette tips**



Example photo of successful run

Tube F Processing

*In **Tube F**, you should have 2 mL from original RPMI + 0.5-1 mL from swab cryopreservative + 1 mL from tube C supernatant = ~4 mL total*

67. Add 1 mL RLT + BME (or Seq-Well complete lysis buffer) to **Tube F**.

68. Distribute contents of **Tube F** into 3 cryovials per sample

You can save time here by labeling cryovials/plates for viral and bulk lysates the day before!

69. Snap freeze on dry ice

70. Store at -80°C

Lysates for Bulk RNA-seq

71. For each lysate, add "Volume Necessary for Each Lysate" from Table 7 to one well of a 96 well PCR plate according to the Lysis Storage Plate map.

Typically, we aim to store 3 lysates per sample. To do so we label 3 plates and fill 1 well per sample in each plate, if cell numbers allow.

72. Seal plates and centrifuge at 400xg for 5 minutes

73. Aspirate media

It may be difficult to see a pellet here, so it's okay to leave some residual volume before resuspending in lysis buffer in next step.

74. Resuspend cells in 50 µL Lysis Buffer (RLT + 1% BME). Mix with pipette. Bubbles are okay here.

75. Seal plate with foil seal and spin down briefly.

76. Place plate on dry ice for 15-20 minutes to snap-freeze lysate.

78. Store at -80°C until ready to use.